


AginWhere™ FT



INCI: Bifida Ferment Filtrate

CAS No.: 96507-89-0

EINECS No.: 306-168-4

Technical Data

Item	Standard
Appearance	Transparent liquid
Color	Colorless
Odor	Basically odorless
pH(20°C)	5.0-7.0
Hg	≤ 1ppm
As	≤ 10ppm
Heavy metal(Pb)	≤ 40ppm
Bacterial count	≤ 500cfu/g
Mold and Yeast	≤ 100cfu/g
E.Coli	Negative
Pseudomonas aeruginosa	Negative
Staphylococcus aureus	Negative

Recommendations:

- AginWhere™ FT is clear liquid that can be applied to transparent systems.
- Addition below 50°C, it can be added to water phase to participate in cold emulsification; Thickener is recommended to use.
- Should not be directly added to strong acid and alkali systems.

Features:

- Fermentation, probiotic Bifidobacterium
- Contains metabolites, amino acids, metabolic acids, vitamins
- Increase skin resistance to a variety of external factors

Packaging&storage:

- Generally, AginWhere™ FT is packed in 1kg or 5kg fluorinated bottles
- AginWhere™ FT is recommended to stored in dry and cool place(4~10°C). Strong light and heat should be avoid as well.
- AginWhere™ FT has 1 year shelf life.



No.35 Dingxing Road, Food Industrial Zone, Yangzhou, Jiangsu, 225004, P.R. China

Tel: 0514-8085 2111 Fax: 0514-8581 5959

Email: fci@joyvofci.com

The information in this document is believed to be reliable. Customers are urged to make their own tests to determine if the product meets your specific needs. Nothing contained herein is intended as a guarantee to use our product. In no events the described materials, information, data in the document shall be liable as a part of sale conditions.